# **Sampling of Spider Mites**

## **Protocol 1 – Infested Leaves**

#### **Objective:**

To collect and prepare leaf samples visibly infested with spider mites for later examination or laboratory processing.

#### **Materials**

- Clean paper envelopes or plastic bags (zip-lock)
- Hand lens (optional, for easier viewing)
- Permanent marker and labels
- Cooler box with ice packs (for transport, optional)

#### **Procedure**

- 1. **Select plants:** Identify plants showing typical spider mite symptoms speckling, chlorosis, bronzing, or webbing on leaves.
- 2. **Choose leaves:** Select 3–5 leaves per plant, preferably from different positions (top, middle, lower canopy).
- 3. Collect leaves: Gently pull leaves with your fingertips to avoid dislodging mites.
- 4. **Inspect briefly:** Use a hand lens, if available, to confirm the presence of mites.
- 5. Package samples:
  - Place each sample (per plant or per site) in a labeled paper envelope or plastic bag.
  - o Label with sample ID, date, location, host plant species, and collector name.

## 6. Storage/Transport:

- o Keep samples cool and shaded during transport (avoid direct sunlight).
- o If processing will happen within 24 hours, store the sample at 4 °C or include ice packs in the shipping box; otherwise, continue with preservation (Protocol 2).

# **Protocol 2 – Preservation in Ethanol**

## **Objective:**

To preserve spider mite specimens collected from leaves for molecular analysis.

#### **Materials**

- Sterile tubes of 1.5–2 mL (preferred), or small, capped containers up to 50 mL.
- 70–100% grade ethanol
- Fine brush
- Labeling tape or waterproof marker
- Optional: dissecting microscope

#### **Procedure**

#### 1. Transfer mites:

- o Under a stereomicroscope or with a hand lens, gently remove mites from leaves using a fine brush.
- o Place them directly into a labeled tube containing ethanol.

### 2. Label samples clearly:

o Include sample ID, date, host plant, and collection site.

#### 3. Storage:

o Store tubes at −20 °C for long-term preservation (preferred) or 4 °C for short-term storage.